

E-Z Split Key[®] Cup

Urine Drug Test



Flat Panel for
Photocopying



Instant Results



Key Activation



With its unique dual chamber design, E-Z Split Key Cup urine drug test kit is a benchmark for reliable drug screening. If a sample needs to be preserved before testing, E-Z Split Key Cup will not activate until the integrated key is inserted. The cup's design also showcases a flat side for ease of photocopying and scanning. This innovative feature makes confirmation testing and reporting a breeze. E-Z Split Key Cup can screen for up to 14 common drugs of abuse, and is offered in 17 different configurations to suit your business needs.

TESTS FOR:

Amphetamine (AMP), Barbiturates (BAR), Benzodiazepines (BZO), Buprenorphine (BUP), Cocaine (COC), Ecstasy/Molly (MDMA), Marijuana (THC), Methadone (MTD), Methamphetamine (mAMP), Morphine/Opiates (OPI), Oxycodone (OXY), Phencyclidine (PCP), Propoxyphene (PPX), Tricyclic Antidepressants (TCA)

PRODUCT FEATURES:

- Simple Twist of a Key Initiates Test
- Flat Panel for Ease of Reporting
- SAMHSA Compliant, CLIA-Waived, FDA 510(k) Approved
- Detects Metabolites for Up to 14 Common Drugs of Abuse, Plus Adulterants
- 99% Accuracy, Results Within 5 Minutes
- 25 Cups Per Box

WHY E-Z SPLIT KEY CUP?

E-Z Split Key Cup allows collectors to withdraw samples and initiate testing whenever convenient with a simple twist of a key. Reporting of results is rapid and easy with E-Z Split Key Cup's flat panel. These unique features speed up urinalysis, creating a hassle-free experience for donors and collectors alike. For reliable results that suit your schedule, try E-Z Split Key Cup.

INSTRUCTIONS

1. Tear the foil bag open, remove test cup. Label the device with donor information.
2. Open test cup lid. Urinate directly into the test cup. Be sure to fill the test cup with urine sample between a minimum of 30 ml to a maximum of 110 ml. The minimum and maximum levels are marked on the cup.
3. After urine specimen has been collected, close the lid securely and return cup to specimen collector.
4. For specimen collector: Peel off label to reveal test result. Read results at 5 minutes. Results should not be read after 10 minutes.

INTERPRETATION OF RESULTS

NEGATIVE:

Two color lines appear. One color line in the Control (C) region and another line next to the test (T) region.



POSITIVE:

One color line appears in the control (C) region and no line appears in the test (T) region.



INVALID:

Control line fails to appear. In case of invalid results, review procedure and repeat test using a new test device. If problem persists, discontinue using the lot immediately and contact your supplier.



IMPORTANT INFORMATION

Do not use after the expiration date.

The test device should remain in the sealed pouch until use.

Results are preliminary and positive results do not always mean that drugs are present in the donor. Certain foods and medicines, diet plan drugs, and nutritional supplements may yield a false positive test result with this product. Similarly, if the sample is diluted or the sample is adulterated, the product may yield negative results and not detect a particular drug. Laboratory tests confirm the preliminary results.



TEST	DESCRIPTION	CUT-OFF LEVEL
AMP	AMPHETAMINE	1000 nG/mL
BAR	BARBITURATES	300 nG/mL
BZO	BENZODIAZEPINES	300 nG/mL
BUP	BUPRENORPHINE	10 nG/mL
COC	COCAINE	300 nG/mL
THC	MARIJUANA	50 nG/mL
MTD	METHADONE	300 nG/mL
mAMP	METHAMPHETAMINE	1000 nG/mL
MDMA	ECSTASY/MOLLY	500 nG/mL
MOR/OPI	MORPHINE/OPIATES	300/2000 nG/mL
OXY	OXYCODINE	100 nG/mL
PCP	PHENCYCLIDINE	25 nG/mL
PPX	PROPOXYPHENE	300 nG/mL
TCA	TRICYCLIC ANTIDEPRESSANTS	1000 nG/mL