

# iScreen OFD<sup>®</sup> Oral Fluid Test Kit

Instant Oral Drug Testing Results



Cotinine Testing  
Available



Non-Invasive



Difficult to  
Adulterate



For a comprehensive saliva drug test that yields lab-grade results, iScreen is the perfect solution. Unlike other drug tests, iScreen is a non-invasive test, with a significantly lower risk of adulteration than traditional urine testing. iScreen can identify up to 6 common drugs of abuse with a separate option to test for cotinine, a nicotine metabolite. Once the swab is adequately soaked, results are apparent within 10 minutes. iScreen is an innovative way to expedite workplace drug testing.

## TESTS FOR:

Amphetamine (AMP), Cocaine (COC), Nicotine/Cotinine (COT), Marijuana (THC), Methamphetamine (mAMP), Opiates (OPI), Phencyclidine (PCP)

## PRODUCT FEATURES:

- Tests for Up to 6 Common Drugs of Abuse, or Cotinine (Nicotine Metabolite)
- 1, 5 or 6 Panel Configurations
- Ability to Detect Drug Use Hours Before Urine Testing
- Non-Invasive Testing and No Restroom Required
- 99% Accuracy, Results Within 10 Minutes
- 25 Tests Per Box



2. Bring the pouch to room temperature before opening. Remove the device from the sealed pouch and use the test immediately.
3. Remove cap from the collection sponge. Insert sponge into mouth.
4. Actively swab the inside of the mouth and tongue with the sponge. As soon as the sponge softens slightly, gently press the sponge between the tongue and teeth to ensure complete saturation. The sponge is saturated when no hard spots are felt. Collect sample for at least 3 minutes before removing sponge.

#### COLLECTOR:

1. Gently plunge the OFD sponge into the cap 3 times.
2. Align the arrow on the device with either one of the marks on the cap. Insert the collector vertically into the cap and press down firmly until the cap reaches the thread.
3. Twist the handle clockwise 180 degrees to tighten the cap until the arrow lines up with the other white mark.
4. Place the device horizontally on a clean and level surface with results window facing up. The test will begin. Within seconds you will notice a pinkish migration or flow of sample as the test strips are activated.
5. Wait for colored lines to appear. Read test results at 10 minutes. Do not interpret results after 1 hour to avoid false results.

#### INTERPRETATION OF RESULTS

##### NEGATIVE:

Two color lines appear. One color line in the Control (C) region and another line next to the test (T) region.



##### POSITIVE:

One color line appears in the control (C) region and no line appears in the test (T) region.



##### INVALID:

Control line fails to appear. In case of invalid results, review procedure and repeat test using a new test device. If problem persists, discontinue using the lot immediately and contact your supplier.



#### IMPORTANT INFORMATION

Do not use after the expiration date.

The test device should remain in the sealed pouch until use.

Results are preliminary and positive results do not always mean that drugs are present in the donor. Certain foods and medicines, diet plan drugs, and nutritional supplements may yield a false positive test result with this product. Similarly, if the sample

| TEST | DESCRIPTION         | CUT-OFF LEVEL |
|------|---------------------|---------------|
| AMP  | AMPHETAMINE         | 50 NG/ML      |
| COC  | COCAINE             | 20 NG/ML      |
| COT  | COTININE (NICOTINE) | 30 NG/ML      |
| mAMP | METHAMPHETAMINE     | 50 NG/ML      |
| OPI  | OPIATES             | 40 NG/ML      |
| PCP  | PHENCYCLIDINE       | 10 NG/ML      |
| THC  | MARIJUANA           | 12 NG/ML      |

#### WHY ISCREEN?

iScreen is the perfect solution to quickly test for the most common drugs of abuse. The simplicity of the saliva swab permits testing from any location - without the need for a restroom. The iScreen OFD cotinine test is great for organizations with a zero tolerance smoking policy, or for insurance tobacco testing. Staffing, recruiting and workplace firms can adopt iScreen for pre-employment drug screening needs.

#### INSTRUCTIONS

##### DONOR:

1. Prior to testing, DO NOT place anything in the mouth for at least 10 minutes, including food, drink, gum, tobacco products or other substances.