



URINE



SALIVA



HAIR

Physicians | Hospitals | Government | Employers | Staffing | Addiction Treatment Centers | Shelters | Resellers

## RAPID DRUG TEST CUT-OFF & DETECTION TIME

### MARIJUANA

Cannabis, Weed, Grass, Mary Jane

50  
ng/mL

3-10  
days



12  
ng/mL

up to  
14  
hrs



1  
pg/mg

up to  
90  
days



300  
ng/mL

1-2  
days



20  
ng/mL

up to  
24  
hrs



500  
pg/mg

up to  
90  
days



### COCAINE

Powder, Coke, Dust, Blow, Bump

300  
ng/mL

1-3  
days



40  
ng/mL

up to  
72  
hrs



200  
pg/mg

up to  
90  
days



### OPIATES

Morphine, Avinza®, Duramorph®, Embeda®  
Kadian®, MS Contin®, MSIR®, Roxanol®

1000  
ng/mL

2-4  
days



50  
ng/mL

up to  
72  
hrs



500  
pg/mg

up to  
90  
days



### AMPHETAMINE

Adderall®, Benzedrine®, Dexedrine®, DextroStat®

### METHAMPHETAMINE

Ampedroxyn®, Desfedrin®, Desoxyn®, Methadrine®

1000  
ng/mL

3-5  
days



50  
ng/mL

up to  
72  
hrs



500  
pg/mg

up to  
90  
days



100  
ng/mL

1-3  
days



50  
ng/mL

up to  
48  
hrs



200  
pg/mg

up to  
90  
days



### OXYCODONE

Endocet®, OxyContin®, OxyIR®, Percocet®  
Percodan®, Roxicet®, Roxicodone, Tylox®

# RAPID DRUG TEST CUT-OFF & DETECTION TIME

**MDMA**  
Ecstasy, Molly

500 ng/mL

1-3 days



50 ng/mL

up to 48 hrs



500 pg/mg

up to 90 days



300 ng/mL

3-7 days



50 ng/mL

up to 48 hrs



100 pg/mg

up to 90 days



**BENZODIAZEPINES**  
Ativan®, Emotavil®, Klonopin®, Lorax®, Niravam®, Serax®, Valium®, Xanax®, Xanor®

**BUPRENORPHINE**  
Buprenex®, Butrans®, Suboxone®, Subutex®

10 ng/mL

3-6 days



10 ng/mL

up to 48 hrs



500 pg/mg

up to 90 days



300 ng/mL

4.5-7 days



300 ng/mL

up to 48 hrs



500 pg/mg

up to 90 days



**BARBITURATES**  
Amytal®, Belladonna®, Butibel®, Butisol®, Fioricet®, Fiorinal®, Luminal®, Nembutal®

**METHADONE**  
Dolophine®, Diskets®, Methodose®


300 ng/mL

1-7 days



75 ng/mL

up to 48 hrs



500 pg/mg

up to 90 days



1000 ng/mL

2-7 days



NA

up to 48 hrs



NA

NA



**TRICYCLIC ANTIDEPRESSANTS**  
Adapin®, Aventyl®, Elavil®, Endep®, Ludiomil®, Pamelor®, Norpramin®, Sinaquan®, Tofranil®

**PHENCYCLIDINE**  
Angel Dust, Sernyl

25 ng/mL

7-14 days



10 ng/mL

up to 72 hrs



300 pg/mg

up to 90 days



NA = Not Available for this Specimen  
These general guidelines give approximate detection periods for each drug substance. Interpretation of retention time must take into account variability of urine specimens, drug metabolism and half-life, metabolic rate, fluid intake, body mass, age, overall health and pH.  
\* Hair testing is laboratory based.