

SalivaConfirm Premium



ITEM	DESCRIPTION
HE-SWI-254	5 Panel SalivaConfirm Premium Drug Test with Saturation Indicator; COC/THC/AMP/mAMP/OPI/PCP
HE-SWI-264	6 Panel SalivaConfirm Premium Drug Test with Saturation Indicator; AMP/COC/mAMP/OPI/THC/BZO
HE-SWI-3964	6 Panel SalivaConfirm Premium Drug Test with K2 with Saturation Indicator; AMP/COC/mAMP/OPI/THC/K2
HE-SWI-17104	10 Panel SalivaConfirm Premium Drug Test with Saturation Indicator; AMP/OPI/THC/COC/PCP/OXY/BZO/BAR/mAMP/BUP
HE-SWI-3124A	12 Panel SalivaConfirm Premium Drug Test with Saturation Indicator; AMP/mAMP/COC/OPI/THC/PCP/BZO/OXY/MTD/BAR/BUP/ALCO
HE-SWI-5124A	12 Panel SalivaConfirm Premium Drug Test with Fentanyl with Saturation Indicator; AMP/BAR/BUP/BZO/COC/mAMP/MTD/OPI/OXY/THC/FEN/ALC
HE-SWI-5144A	14 Panel SalivaConfirm Premium Drug Test with Saturation Indicator; AMP/mAMP/COC/OPI/THC/PCP/BZO/OXY/MTD/BAR/BUP/K2 SPICE/FENTANYL/ALCO

DRUG CUT-OFF LEVELS

TEST	DESCRIPTION	CUT-OFF LEVEL
ALC	ALCOHOL	>0.02% B.A.C.
AMP	AMPHETAMINE	50 NG/ML
BAR	BARBITURATES	300 NG/ML
BZO	BENZODIAZEPINES	50 NG/ML
BUP	BUPRENORPHINE	10 NG/ML
COC	COCAINE	20 NG/ML
FEN	FENTANYL	10 NG/ML
K2/SPICE	K2 SPICE	5 NG/ML
THC	MARIJUANA	11-NOR- $\Delta(9)$ -THC-9 COOH - 12NG/ML $\Delta(9)$ -THC - 75 NG/ML
mAMP	METHAMPHETAMINE	50 NG/ML
MOR/OPI	MORPHINE/OPIATES	40 NG/ML
MTD	METHADONE	75 NG/ML
OXY	OXYCODONE	50 NG/ML
PCP	PHENCYCLIDINE	10 NG/ML