



DrugConfirm™ Advanced

TOP-SELLING
INTEGRATED URINE
DRUG TEST CUPS

DrugConfirm™ Advanced is the only urine drug test cup that can screen up to 16 common drugs at once. This industry-preferred drug screening solution delivers up to 99% accurate results in minutes. Our trademarked DrugConfirm™ Advanced urine test cups also offer built-in adulterants to prevent tampering of urine specimens. See which DrugConfirm™ Advanced configuration is the right fit for your drug testing protocol:



For **Clinical, Government,
and Workplace Professionals**




Prevent **Specimen
Adulteration**



CLIA Waived
Options Available



Low Cut Offs for
Sensitive Detection



PRODUCT FEATURES :

- CLIA Waived, FDA 510(K) Cleared Available
- Screen up to 16 Drugs at Once
- Built-In Temperature and Specimen Validity Strips
- Rapid Results in Minutes
- Option for Custom Configurations
- Up to 99% Accuracy
- For Clinical, Workplace, and Forensic Use
- 25 Cups Per Box

TESTS FOR :

Alcohol (EtG)*, Amphetamine (AMP), Barbiturates (BAR), Benzodiazepines (BZO), Buprenorphine (BUP), Cocaine (COC), Ecstasy/Molly (MDMA), **Fentanyl (FEN)***, Marijuana (THC), Methadone (MTD), Methamphetamine (mAMP), Morphine/Opiates (MOR/OPI), Oxycodone (OXY), Phencyclidine (PCP), **Synthetic Marijuana (K2/Spice)***, **Tramadol (TRA)***, Tricyclic Antidepressants (TCA), Adulterants Available: Creatinine (CR), Specific Gravity (SG), Oxidants (OX)

***For Forensic Use Only**

CB00083118

DRUGCONFIRM™ ADVANCED CUP CONFIGURATION OPTIONS

Drug Cut-Off Concentration Levels	12-Panel CLIA Waived (#6125)	14-Panel CLIA Waived (#1144)	12-Panel 2.0 (#8124)	13-Panel 3.0 (#8134)	16-Panel (#6164KT)
Alcohol (EtG)* 300/500 ng/mL		●	●		●
Amphetamine (AMP) 1000 ng/mL	●	●	●	●	●
Barbiturates (BAR) 300 ng/mL	●			●	●
Benzodiazepines (BZO) 300 ng/mL	●	●	●	●	●
Buprenorphine (BUP) 10 ng/mL	●	●	●	●	●
Cocaine (COC) 300 ng/mL	●	●	●	●	●
Ecstasy/Molly (MDMA) 500 ng/mL	●	●	●	●	●
Fentanyl (FEN)* 200 ng/mL		●	●		●
Marijuana (THC) 50 ng/mL	●	●	●	●	●
Methadone (MTD) 50 ng/mL	●	●	●	●	●
Methamphetamine (mAMP) 1000 ng/mL	●	●	●	●	●
Morphine/Opiates (MOR/OPI) 300/2000 ng/mL	●	●	●	●	●
Oxycodone (OXY) 100 ng/mL	●	●	●	●	●
Phencyclidine (PCP) 25 ng/mL	●			●	●
Propoxyphene (PPX) 300 ng/mL	●			●	
Synthetic Marijuana (K2/Spice)* 50 ng/mL		●	●		●
Tramadol (TRA*) 200 ng/mL					●
Tricyclic Antidepressants (TCA) 1000 ng/mL				●	
Adulterants Available (CR, SG, OX)	●	●	●	●	

CB00083118