

DRUGCONFIRM™ ADVANCED CUP CONFIGURATION OPTIONS

Drug Cut-Off Concentration Levels	12-Panel CLIA Waived (#6125)	14-Panel CLIA Waived (#1144)	12-Panel 2.0 (#8124)	13-Panel 3.0 (#8134)	16-Panel (#6164KT)
Alcohol (EtG)* 300/500 ng/mL		●	●		●
Amphetamine (AMP) 1000 ng/mL	●	●	●	●	●
Barbiturates (BAR) 300 ng/mL	●			●	●
Benzodiazepines (BZO) 300 ng/mL	●	●	●	●	●
Buprenorphine (BUP) 10 ng/mL	●	●	●	●	●
Cocaine (COC) 300 ng/mL	●	●	●	●	●
Ecstasy/Molly (MDMA) 500 ng/mL	●	●	●	●	●
Fentanyl (FEN)* 200 ng/mL		●	●		●
Marijuana (THC) 50 ng/mL	●	●	●	●	●
Methadone (MTD) 50 ng/mL	●	●	●	●	●
Methamphetamine (mAMP) 1000 ng/mL	●	●	●	●	●
Morphine/Opiates (MOR/OPI) 300/2000 ng/mL	●	●	●	●	●
Oxycodone (OXY) 100 ng/mL	●	●	●	●	●
Phencyclidine (PCP) 25 ng/mL	●			●	●
Propoxyphene (PPX) 300 ng/mL	●			●	
Synthetic Marijuana (K2/Spice)* 50 ng/mL		●	●		●
Tramadol (TRA*) 200 ng/mL					●
Tricyclic Antidepressants (TCA) 1000 ng/mL				●	
Adulterants Available (CR, SG, OX)	●	●	●	●	

CB00083118